

Abstract

Protein-based and cellular therapeutics have revolutionized the treatment of cancer and other diseases. However, the need to produce complex biologics with limited resources and shortened timelines has created a demand for technologies that can express multiple therapeutic moieties in a range of cell types and culture processes. This poster describes the use of MaxCyte's flow electroporation-based delivery platform for rapid, transient production of monoclonal mouse and human antibodies, Fc-fusions and other proteins in a range of host cells. Emphasis will be placed on the use of large-scale transient transfection for optimizing media and feed strategies in batch fed cultures using multiple CHO cell lines. Data on protein titers and cell viability will be shared to illustrate how a flexible transfection platform enables the development of optimized production processes.

High Titer, CHO-S IgG Production Using Media of Choice
Media & Feed Optimization Produces 7x Improvement in Antibody Titers





Figure 1: 7x Improvement in Titer Upon Optimization of Media and Feed. 2e9 CHO-S cells were transfected with a human IgG plasmid via a single large-scale electroporation. 1.2e8 cells were seeded in each of 14,125mL flask. 30mL of media (Medium 1, 2, 3, 4 or 5) was added per flask bringing the initial seeding density to 4e6 cells/mL. Cells were either grown as batch or fed-batch cultures following the feed schedule indicated in the table. Antibody titers, viability (data not shown) and metabolic profiles (data not shown) were analyzed on Days 4, 6, 8, 10 and 12 post transfection.

Media Flexibility Enables Significant Cost Savings
1.8x Increase in Fc-fusion Titers With a 75% Decrease in Media Cost




Figure 2: Significant Reduction in Media Costs While Increasing Fc-fusion Production. 2e9 ExpiCHO cells were transfected with an Fc-fusion protein expression plasmid via a single large-scale electroporation. 1.5e8 cells were seeded per 125mL flask with 25mL of one of three commercially available media and one of four feeds. Cells were either grown as batch or fed-batch cultures following the feed schedule indicated in the table. Fc-fusion protein titers were analyzed on Days 9 post transfection.

Media & Feed Flexibility for Glyco-pattern Engineering
Production of Differently Glycosylated IgG from CHOZN Cells




Figure 3: Media and Feed Flexibility Post Electroporation Allow Customization of Therapeutics. CHOZN® GS-/- cells were transfected with an hlgG expression plasmid via 3 independent electroporations using the MaxCyte STX. Cells from independent electroporations were cultured post electroporation in ExCell CD Fusion medium using two different feed supplements. Antibody glycosylation patterns were analyzed from 8 day fed-batch cultures.

Increasing Laboratory Productivity Using Flow Electroporation Technology & Higher Seed Densities

4x to 28x Titer Improvements for Human & Mouse IgG While Cutting Production Time in Half




Figure 4: Increased Productivity Through Seed Density Optimization. 6e9 ExpiCHO cells were transfected with either a human or mouse IgG plasmid via a single large-scale electroporation. Transfected cells were seeded in 125mL flasks at three different densities: low (4e6 cells/mL), medium (8e6 cells/mL) or high (1.2e7 cells/mL). Cells were fed as indicated in the table. Antibody titers and cell viability were analyzed on Days 4, 5, 6 and 7 post transfection. Results were compared to the customer's established ExpiCHO process using ExpiCHO medium.

Seed Density Optimization for Improved Antigen Titers

2.5x Higher Titers Using Flow Electroporation Technology & Increased Cell Densities




Figure 5: Higher Antigen Expression Levels Upon Switching to MaxCyte Electroporation & Optimizing Seed Density. 2e9 ExpiCHO cells were transfected with an antigen expression plasmid via a single large-scale electroporation. Transfected cells were seeded in CD OptiCHO medium at three different densities: low (6e6 cells/mL), medium (1.0e7 cells/mL) or high (1.4e7 cells/mL). Cells were fed on Day 5 and Day 9 with commercially available supplements. Antibody titers and cell viability were analyzed at various times over the course of a 14 day culture. Results were compared to the customer's established ExpiCHO process using ExpiCHO medium.

Increased Productivity Using Lower Cost Media & Feed
IgG Titers Improved to >4g/L




Figure 6: Higher IgG Titers Using Flow Electroporation Technology & Lower Cost Media and Feed. ExpiCHO cells adapted to CDM4MAB media prior to transfection. 2e9 cells were transfected with a IgG expression plasmid using a single large-scale electroporation. Transfected cells were seeded at either 6e6 cells/mL or 1e7 cells/mL in CDM4MAB media. Cultures were supplemented every other day with specified feeds and additional glutamine. Antibody titers were determined at various times throughout a 13 day culture. Two independent ExpiCHO transfections were performed using vendor reagents and cultured and fed according to vendor's 'Max Titer' protocol.

Summary

- MaxCyte's Flow Electroporation™ Technology enables high efficiency, high viability transfection of a variety of CHO cell lines, including ExpiCHO cells.
- MaxCyte scalability allows the use of a single large-scale electroporation for in-depth process development, eliminating multiple transfection as a study variable.
- The MaxCyte platform allows the use of any media or culture supplements enabling:
 - significantly lower consumable costs (75% reduction in price enabled) with equal to or improved protein titers in the cell line of choice
 - use of specialty media such as those used to direct glycosylation patterns
- MaxCyte transfection enables seed densities >1e7 cells/mL post electroporation while maintaining high antibody expression, thereby improving laboratory productivity.
- Process optimization enabled substantially higher expression levels which allows shorter production runs.
- High viability and transfection efficiency results in strong expression of more difficult-to-express proteins such as Fc-fusion proteins and mouse IgGs.