



# Enabling Scalable Cell Engineering Using the MaxCyte<sup>®</sup> Electroporation Platform


Andrew Mancini, PhD  
Field Applications Scientist  
March 28th, 2022

 MaxCyte<sup>®</sup> e∞per<sup>†</sup> AT<sub>x</sub> GT<sub>x</sub> ST<sub>x</sub>

are registered trademarks of MaxCyte, Inc.




*This document may contain proprietary content and/or confidential information owned by MaxCyte, Inc. and/or third parties that have granted MaxCyte, Inc. limited rights to use their content/information. Your receipt of this document does not constitute a license or right to use, copy, distribute, disclose, derive, or otherwise use any of the content/information contained herein. MaxCyte makes no warranties or representations regarding the accuracy, authenticity, or origin of any of the content/information in this document except as expressly stated herein.*

# High Performance Electroporation Technology *from Concept to Clinic*



MaxCyte enables significant reduction in **Timelines, Risk** and overall **Cost**

# MaxCyte Electroporation Instruments

Instrument	Touch Screen	LED Indicators	Static EP	Flow EP	Barcode Reader	21CFR	Master File	cGMP
 <span data-bbox="351 464 624 571">GTx<sup>®</sup></span>	●	●	●	●	●	●	●	●
 <span data-bbox="364 742 624 856">STx<sup>®</sup></span>	●	●	●	●				
 <span data-bbox="377 1021 624 1135">ATx<sup>®</sup></span>	●	●	●					

# The Leading Choice For Advanced Cell Engineering



## High-performance engine accelerating the clinical translation of next generation cell therapies

- Twenty years of experience in primary cell transfection, protein production, and genome engineering
- Continued expansion of cell therapy partnerships with leading industry innovators and [multiple programs currently in clinical trials](#)
- Trusted by the top academic, translational and commercial groups for their R&D, PD and MFG work



# High Transfection Efficiency with Any Cargo



Plasmid DNA



Linear dsDNA




ssODN



mRNA



Antisense Oligonucleotides (ASOs)




Small RNAs




CRISPR RNP




Recombinant Protein



Cell Lysates



# Express Any Molecule with the MaxCyte Expert™ Platform



Ion channels  
GPCRs  
Transporters  
Receptors




CRISPR-Cas9  
ZFNs  
TALENs  
siRNAs/ASOs  
Base and Prime Editors



Antibodies  
Bi-specifics  
Fusion Proteins  
Antigens



Transcription  
Factors  
  
Epigenetic  
Regulators



Enzymes  
Reporters



Lentivirus  
Alphavirus  
Adenovirus/AAV  
VLPs  
Subunit vaccines




## Immortalized Mammalian Cells

10T½	L6
C127	LLC-PK1
CAP/CAP-T	MDCK
Ba/F3	NIH/3T3
BHK-21	NRK
C2C12	PC12
CHO	Per.C6
COS-7	RAW264.7
HEK-293	RBL
Jurkat	Vero

## Mammalian Cancer Cell Lines

A549	Neuro2a
B16	NS0
EL4	PANC-1
HCT-116	Renca
HeLa	SK-N-SH
HT-1080	SK-OV-3
K562	SNU-1
LNCaP	Sp2/0
MCF7	THP-1
MG-63	YB2/0



## Hematopoietic Cells

HSCs	Dendritic cells
PBMCs	NK cells
B cells	Granulocytes
T cells	Monocytes

## Stem Cells

iPSCs
Embryonic stem cells
Mesenchymal stem cells
Neural stem cells
Neural progenitor cells

## Primary Cells

Myoblasts
Keratinocytes
Neurons
Fibroblasts
Cardiomyocytes

## Insect Cell Lines

S2
SL3
Sf9
Sf21

**+ 10 optimization protocols for new cell types and applications**



# Seamless Scalability Using MaxCyte Processing Assemblies




Feature	OC-25x3	R-50x3	R-50x8	OC-100x2	OC-100	OC-400	R-/G-1000	CL-1.1	CL-2
PA Type									
High Vol.	25 $\mu$ L	55 $\mu$ L	55 $\mu$ L	100 $\mu$ L	100 $\mu$ L	400 $\mu$ L	1 mL	3.5 mL	100 mL
Low Vol.	15 $\mu$ L	45 $\mu$ L	45 $\mu$ L	50 $\mu$ L	50 $\mu$ L	200 $\mu$ L	400 $\mu$ L	1 mL	10 mL
# Samples	3	3	8	2	1	1	1	1	1
High Cell	$5 \times 10^6$	$1 \times 10^7$	$1 \times 10^7$	$2 \times 10^7$	$2 \times 10^7$	$8 \times 10^7$	$2 \times 10^8$	$7 \times 10^8$	$2 \times 10^{10}$
Low Cell	$7.5 \times 10^4$	$2.25 \times 10^5$	$2.25 \times 10^5$	$2.5 \times 10^5$	$2.5 \times 10^5$	$1 \times 10^6$	$2 \times 10^6$	$5 \times 10^6$	$5 \times 10^7$
ATx	➔								
GTx	➔								
STx	➔								





# A Fully Optimized CRISPR Workflow for Drug Discovery in T Cells

SLICE-Enabled Whole Genome Screens to Investigate Immune Dysregulation




# Highly Efficient Homozygous Correction of the DYSF Gene in Miyoshi Myopathy Patient iPSCs by ssODN-Mediated Knockin

>70% Homozygous corrected iPSCs enable optimized selection of isogenic clones for disease modeling and drug screening




# A GMP Compatible, Non-Viral CAR T Cell Manufacturing Process






# Highly Efficient CAR Knock-In and Robust T Cell Expansion


40-62% target knock-in



30-fold CAR+ T cell expansion over a 10 day manufacturing process



Robust cancer cell killing





Have questions? Want to learn more?  
Visit Booth #123!

Thank you!

 MaxCyte® e∞pert<sup>™</sup> AT<sub>x</sub><sup>™</sup> GT<sub>x</sub><sup>™</sup> ST<sub>x</sub><sup>™</sup>

are registered trademarks of MaxCyte, Inc. in the U.S.A.